

VENOSAN®

SHEER

Compression Stockings

SOFT & ELEGANT



VENOSAN®

SHEER

Ccl. II
23-32 mm Hg
(3,1-4,3 kPa)

Compression Stockings

Medias de compresión | Compression Stockings | Bas compressifs | Calze compressive




AD
Debajo de la rodilla
Below knee
Bas au mollet
Gambaletti



AGH
Hasta el muslo
Thigh length
Bas à cuisse
Calza alla coscia


Compresión gradual | Graduated compression | Compression dégressive | Compressione graduata



Banda superior de encaje | Lace top band | Adhésive ornée | Bordo autoreggente pizzo




Exterior
Outside
Face externe
Lato esterno



Interior
Inside
Côte intérieur
Lato interno

 **SWISS COMFORT**



Effective Compression

Mejora el flujo sanguíneo
Improved blood flow
Amélioration de la circulation sanguine
Miglioramento del flusso sanguigno

Leg Contouring


Alta eficacia de compresión
Beautiful legs through shaping compression
Compression modelante pour les jambes
Compressione modellante per le gambe

Anatomically Knitted

Para máxima comodidad, no se resbala
For maximum comfort, does not slip
Pour un confort optimal, ne glisse pas
Per il massimo comfort, non scivola

Medidas de circunferencia | Circumference measurements |
 Mesures de circonférence | Misure circonferenza (cm)

	S	M	L	XL
cB	20-23	23-26	26-29	29-31
cC	29-39	33-42	36-45	39-47
cD	28-38	32-41	35-44	38-46
cG	48-63	54-67	58-72	63-74



Longitud | Length | Longueur | Lunghezza (cm)

AD	ID	35-42
AGH	IG	63-78

Artículo | Article | Article | Articolo

TIPO	*COLORES	PUNTERA	S	M	L	XL
AD	Beige	cerrada	SG22011	SG22012	SG22013	SG22014
	Black	cerrada	SN22011	SN22012	SN22013	SN22014

TIPO	*COLORES	PUNTERA	S	M	L	XL
AGH	Beige	cerrada	SG22601	SG22602	SG22603	SG22604
	Black	cerrada	SN22601	SN22602	SN22603	SN22604

*Colores



Beige



Black



Composición | Composition | Composition | Composizione

65% Polyamide | 35 % Elastane





VENOSAN®
Compression Stockings



VENOSAN® is a registered
trademark of SWISSLASTIC AG
ST. GALLEN, Switzerland © 2017
www.venosan.com

SWISSLASTIC AG ST. GALLEN
Rorschacher Strasse 304
CH-9016 St. Gallen
Switzerland