

# Silicoplant

Ortesis de silicona viscoelástica

Viscoelastic silicone orthoses  
Orthese en silicone viscoélastique  
Orthese aus viskoelastischem silikon  
Ortesi in silicone viscoelastico  
Órtese de silicone viscoelástica



**EMO**  
especialistas  
en ortopedia

## Separador de carrete con anillo

Reel toe separator  
with ring

POD201



MADE IN SPAIN



ESPECIALIDADES MÉDICO ORTOPÉDICAS, S.L.  
Pol. Ind. La Pascueta, Camino de La Pascueta, 40  
46200 Paiporta (Valencia) España  
[www.emo.es](http://www.emo.es) +34 963 132 051

TALLAS  
SIZES  
S, M, L

**Bienestar  
para tus pies**  
Alivio y confort.  
Wellness for your feet  
Relief and comfort.



especialistas  
en ortopedia

**Silicoplant**  
Ortesis de silicona  
viscoelástica  
Viscoelastic silicone orthoses

# Silicoplant

Ortesis de silicona viscoelástica

Viscoelastic silicone orthoses  
Orthese en silicone viscoélastique  
Orthese aus viskoelastischem silikon  
Ortesi in silicone viscoelastico  
Ortese de silicone viscoelástica



**EMO**  
especialista  
en ortopedia

**POD201**

## Silicoplant Separador de carrete con anillo

### INDICACIONES

- Hallux valgus flexible.
- Post-cirugía del juanete.
- Callosidades interdigitales ( ojos de gallo o helomas )

### CARACTERÍSTICAS

- Separador del primer dedo en hallux valgus.
- Fabricado con silicona médica extrablanda.
- Se adapta a la temperatura corporal.
- Lavable.

### EFFECTOS

- Alineación pasiva de la primera articulación metatarsofalángica en hallux valgus.
- Flexible para calzados anchos.
- Evita la formación de callosidades interdigitales.
- Disminuye las presiones en el espacio interdigital.

## Silicoplant Reel toe separator with ring

### INDICATIONS

- Hallux valgus Flexible.
- Post-surgery Bunion.

### FEATURES

- Separator toe for hallux valgus.
- Made of extra-soft medical silicone.
- It adapts to body temperature.
- Washable.

### EFFECTS.

- Passive alignment of the first metatarsophalangeal joint in hallux valgus.
- Flexibility for wide shoes.

TODAS LAS ORTESIS DE PODOLOGÍA EMO EN:

ALL PODIATRY ORTHOSES EMO IN:

[www.emo.es](http://www.emo.es) +34 963 132 051

